Artica • Astra • Unika Plus • Dea • T Paper • Traccia • Test • Sterildocs • Autoclavi

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The sterilization cycle

Infection prevention consists of many aspects and its realization is based mainly on knowing the objective to be reached and on the choice and optimum use of the available means. One of these aspects regards sterilization processes. An "absolute" meaning is, by definition, attributed to the final result of a sterilization process but, in fact, it is the realization of a "probability" that closely depends on the methods used to prepare and pack the material, the efficiency of the means and their correct use, the suitability of the measures chosen for preserving and protecting the state that the materials have reached at the end of the process. It follows that, for a correct approach to infection prevention connected with the use of medical devices, it is important to have the necessary resources and instruments.

Operators and persons in charge of the sterilization process are required to keep themselves up to date with the standards that govern it in order to set operating protocols that are in line with the scientific, technological and procedural state of the art and also to be informed about which instruments are necessary for guaranteeing the final objective, which is to produce sterile material.

The sterilization process is a precise set of procedures that must be followed scrupulously and carefully. The incorrect application of the procedures, or their partial adoption, may compromise the entire process and the targeted result. The preparation of medical devices includes cleaning, sub-divided into several stages: decontamination, cleaning and rinsing, followed by drying and packaging of the load.

Tecno-Gaz is the only company able to offer you a complete, rational, ergonomic, safe sterilization programme.

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Aim of the process

This stage is extremely important since it complies with the Italian legislative decree D.Lgs 626/94 on the protection of workers. Decontamination prevents potentially contaminating instruments being handled by the operator. This operation precedes the actual cleaning of the device and its purpose is to remove most of the organic material present on its surface. To decontaminate the devices, they are immersed in a solution containing chemical agents (usually a disinfectant) so that the organic material, with its possible microbial load, is dissolved without the operator having to handle it directly. So that the decontamination procedure is efficacious, the more complex instruments are dismantled or opened, as far as possible, before being immersed, ensuring that the hollow parts are accessible. After decontamination, the medical devices must be rinsed under cold running water.

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Operating Protocols Video

Tecno-Gaz has produced a teaching kit on sterilization procedures for dental assistants. The kit consists of written procedures and DVD with filmed procedures. (procedures-check tests-useful advice)

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Artica vessels Chemical disinfection

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Artica vessels are made for all disinfection liquids. Tecno-Gaz offers you Clorex 4, a concentrated disinfectant detergent solution suitable for instruments and materials. Ideal for use in decontamination vessels and ultrasonic tanks.

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Disinfection vessels are not technological products and it would therefore have been very simple to make the same old device, understood as a container for liquids. But Tecno-Gaz decided to start a real research and development project even on these devices and the result was a truly positive one, since the Artica decontamination and disinfection vessels have unique characteristics and features that guarantee the operator simple, rational and above all safe use.

Artica vessels are made with a special plastic that is compatible with chemical liquids. They have a shaped bottom that allows the instrument container to remain suspended, so preventing the instruments from coming into contact with the residues resting on the bottom of the vessel.

Artica vessels are equipped with:

- full removable lid, which allows the liquid to be changed easily and quickly.
- partial lid, which allows instruments to be inserted without having to open the vessel completely with consequent escape of odours and fumes.
- support equipped with mechanical timer, which allows the time the instruments stay in the vessel to be programmed and checked and guarantees that the disinfection process is completed correctly.
- internal basket made of plastic.

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ARTICA vessels are available in 2 versions:

art. 2086-S-3 Artica 1 Vessel

- Standard size
- Complete with support
- Timer
- Internal vessel made of plastic.

art. 2087-S-3 Artica 2 Vessel

- Large size
- Complete with support
- Timer
- Internal vessel made of plastic in the S-2 version.

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Aim of the process

Cleaning of the devices to be sterilized is an essential requisite for the successful outcome of the whole process. The purpose of cleaning is to reduce the amount of microbial contamination by over 90% and to remove the organic material left behind by the preceding procedure. If these remain on the devices they can hinder the action of the sterilizing agent and defeat the whole process. If an instrument is not clean, its sterilization cannot in any way be obtained. But the sterilization process cannot detect if an instrument is "still dirty" and cannot establish that, as such, it is not suitable for being sterilized. Not only that, the checking of the physical parameters of the corresponding cycle can wrongly validate a process that in reality has not been efficacious on the whole load. **It is therefore up to the professionalism of the operators to guarantee correct cleaning of the medical device.** Cleaning makes use of the chemical action of a detergent, preferably one that is enzyme based and at the same time non-corrosive, which must be used strictly at the concentrations and for the times recommended by the supplier.

It is necessary to renew the solution frequently to avoid accumulation of residues that would corrode the devices and reduce the cleaning action.

Ultrasonic tanks are the most suitable devices for being sure of a complete removal of organic material and hard materials. The cleaning procedure must be carried out in a dedicated room or space. The cleaning operators must use personal protective equipment (goggles-visors-gloves-aprons).

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Cleaning is a priority and crucial stage for the successful outcome of the sterilization process. Tecno-Gaz therefore decided to give maximum attention to ultrasonic tanks, since this is an extremely important device.

The ASTRA ultrasonic tanks have been designed and developed in order to be able to supply a completely electronically controlled automatic instrument that is highly efficient and easy to use. All the applied technical solutions have been adopted to make the operators' work simple and rational.

ASTRA vessels are available in 12 different versions; 4 are standard and the others are made on request.

The standard versions are the following:

- **ASTRA 3 D with electronic control (3 litres)**
- **ASTRA 3 S with analogue/mechanical control (3 litres)**
- **ASTRA 6 D with electronic control (6 litres)**
- **ASTRA 6 S with analogue/mechanical control (6 litres)**

The versions made on request are models with the following litre capacities: 9.5 - 14 - 18 - 21 in both electronic and analogue/mechanical version.

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The sterilization process is a precise set of procedures that must be followed scrupulously and carefully. The incorrect application of the procedures, or their partial adoption, may compromise the entire process and the targeted result.

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Concentrated disinfectant and detergent solution for medical devices, particularly suitable for disinfection vessels and ultrasonic tanks.

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Concentrated detergent for general use. 1 kg pack.

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The analogue/mechanical versions are extremely ductile models, with excellent technical features:

- Setting of time with mechanical timer
- 60°C heating system that can be activated through a control
- Latest generation transducers
- Excellent sound proofing

All the vessels are equipped with accessories: mesh basket, 2 bekers, beker holder, lid.

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The electronic versions are the most complete and technological solutions, with various options:

- Electronic setting of the time and temperature
- Integrated, electronically controlled heater
- Automatic degassing
- New generation transducers
- Liquid drain connector
 - Excellent sound proofing

All the vessels are equipped with accessories: mesh basket, 2 bekers, beker holder, lid, liquid drain tube with tap.

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dimensioni vaschette / dimensions of all tanks

| | | |
|--------------------|---------------|---------------|
| Larghezza Widht | 240 mm | 300 mm |
|--------------------|---------------|---------------|

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| | | |
|---------------------|---------------|---------------|
| Profondità Depth | 140 mm | 240 mm |
|---------------------|---------------|---------------|

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Aim of the process

After cleaning has finished, the instruments must be properly rinsed to remove residues of materials and all traces of detergent that could interact with the sterilizing agents.

After rinsing, the devices must obligatorily be dried to prevent corrosion and white marks; also, water residues can compromise the subsequent sterilization process.

To perform this operation, we recommend the use of mechanical means and do not recommend the use of paper cloths or other manual methods since there is a high risk of being cut or pricked, with serious consequences.

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Demonstrating yet again its great attention to growing operating needs, Tecno-Gaz has produced a unique and extremely rational mechanical device for drying the materials to be sterilized.

This device allows the obligatory drying to be performed and to avoid any handling of instruments especially sharp ones.

Unika is simple, rational and very fast to use. It guarantees excellent drying, which prevents any possibility of corrosion and white marks on the instruments.

Unika is an indispensable and very useful instrument that helps in the preparation of instruments and above all safeguards the operators from being pricked or cut. The average time required to obtain a good drying is only about two minutes.

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Aim of the process

The purpose of packaging is to guarantee that, after sterilization, the devices remain sterile and are protected from contamination that could occur during transport and storage.

The load must be packaged in a way that guarantees the removal of air inside the package and allows contact with the sterilizing agent. The materials most commonly used for packaging are: **paper-polypropylene laminate.**

Paper-polypropylene laminating is preferably used for small devices and allows the contents to be identified. The devices used in surgeries are packed in a single bag while those used for invasive procedures and those intended for operating theatres are packed in a double bag to guarantee use with the "no touch" technique. The paper-polypropylene bags are subjected to positive pressure, and therefore the seal must be at least 10 mm wide so that they do not open during the sterilization cycle.

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Sealing machine mod. Dea

DEA is the new electronically controlled thermosealing machine that allows all paper/plastic laminates to be heat sealed.

Sealing takes place through the actuation of the special ergonomic lever incorporated in the body of the equipment, so that **it takes up as little room as possible.** The shape is modern and pleasing and the structure has been made with the aim of using the roll holder vertically, so that **it takes up as little room as possible and so that it can fit into any cupboard compartment.**

The heating element used for sealing works by direct heating and this **avoids frequent periodic maintenance so saving time and money.**

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Tecno-Gaz supplies not merely equipment, but a sterilization program, a rational and ergonomic concept that helps you to complete all the phases of the sterilization cycle simply and safely.

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Sealing machine mod. Opera Automatic

Opera Automatic is a roll sealing machine, with automatic feed system with photocell. This is the maximum expression of technology applied to medical sealing, and is used in medical centres, dentists' surgeries, hospitals and clinics, ensuring very high output, maximum operating speed and no routine maintenance, and uses precut bags sealed on 3 sides (these are available in various sizes and can be supplied by Tecno-Gaz).

Opera Automatic is complete with integrated printer that, during sealing, automatically stamps the packaging date and expiry date, as well as the required regulatory symbols, on the bag. This very important solution allows you to always have control over the real sterility maintenance state of the instruments. On the side of the equipment there is also an electronic control and management card that allows the parameters of use to be adjusted. The use of this machine is simple and rational. All the operator has to do is insert the bag with the instrument on the sealing machine support and the photocell automatically and autonomously feeds the bag, seals it and stamps the packaging data. This action can be continuative and repetitive to ensure very high output and considerable reduction of times. All the practices that plan to certify and document all the sterilization procedures must equip themselves with **Opera Automatic**.

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Tracking

The tracking of the instruments is a basic procedure for the handling of sterilized materials. Tracking guarantees the recording of the “path” taken by each and every instrument, from preparation to sterilization, until it is used. This is made possible thanks to a data cross checking management system. Tracking is an activity that forms an integral part of the sterilization cycle procedure and should be implemented after the wrapping phase and before the autoclave is loaded. This procedure is indispensable in order to guarantee maximum legal protection to the operators and to the health director.

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
The “Traccia” system

There are various instrument tracking systems, all of which are managed by technologically advanced information processing software, bar code readers, optic pens etc... These systems are valid and efficacious, but only suitable for computerized structures. In the case of environments that are not totally computerized, objective difficulties frequently exist that enter into conflict with the organization of the surgery. This is why Tecno-Gaz has designed a new and exclusive system called “TRACCIA”, specifically researched and developed for use in a surgery context, ideal for all types of facilities. TRACCIA has been designed so as not to generate problems and obstacles for operators in the overall management of instrument and sterilized materials tracking. It is an extremely valid, simple and practical system which guarantees absolute tracking of sterilized materials.


The Traccia guarantees:

- Tracking of the sterilized load
- Tracking of each individual instrument and identification of the patient on which it was used
- Management of the instruments and control of the sterility maintenance period.

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
Adhesive label of the Traccia system with the information necessary for the tracking of the sterilization.

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The labels of Traccia system change colour after sterilization to grant the result of the cycle.

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T Paper

Tecno-Gaz spa offers a range of sterilization bags and rolls called T paper.

These devices are used in the medical, hospital, dentistry, veterinary, beauty and podology fields and by all operators who use steam or gas sterilization.

The **T paper** products are made with plastic laminate (polyester-polypropylene) and special medical paper, which is moisture resistant with high protection from bacteria.

The sealing is multiline, carried out hot, which prevents air bubble formation, ensures perfect closing and does not allow dust particles to enter the bag, so guaranteeing perfect maintenance of sterility.

The inks used are water based, and there are colour change indicators every 10 cm on the rolls and individually on each bag.

The indicators are for colour change with steam or gas.

All the devices are manufactured according to directive EN 868.

Product info

| | CARTA STANDARD/STANDARD PAPER | CARTA SPECIALE/SPECIAL PAPER |
|-------------------------------------------------------------|---------------------------------------------------------------------|---------------------------------------------------------------------|
| Lunghezza rotoli Roll length | 200 m (standard) 100 m (a richiesta/on demand) | 200 m (standard) 100 m (a richiesta/on demand) |
| Colore del materiale plastico Colour of plastic material | verde standard/standard green azzurro a richiesta/blue on demand | verde standard/standard green azzurro a richiesta/blue on demand |
| gr m ² | 60 gr m ² | 70 gr m ² (a richiesta/on demand) |


T Paper
buste autosigillanti
self-sealing bags

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| CODICE CODE | ROTOLI STERILIZZAZIONE 200 m STERILIZATION PAPER ROLLS 200 m | | |
|-------------|-------------------------------------------------------------------------------------------------------|-------|-----------------------------------------------------------|
| | | 516/S | 100X200mm conf. 1000 pz. 100X200mm 1000 pieces package |
| 400/S | rotoli 50 mm x 200 m conf. 8 pz. rolls 50 mm x 200 m 8 piece package | 518/S | 100X250mm conf. 1000 pz. 100X250mm 1000 pieces package |
| 402/S | rotoli 75 mm x 200 m conf. 8 pz. rolls 75 mm x 200 m 8 piece package | 520/S | 100X300mm conf. 1000 pz. 100X300mm 1000 pieces package |
| 404/S | rotoli 100 mm x 200 m conf. 4 pz. rolls 100 mm x 200 m 4 piece package | 532/S | 150X200mm conf. 1000 pz. 150X200mm 1000 pieces package |
| 406/S | rotoli 150 mm x 200 m conf. 4 pz. rolls 150 mm x 200 m 4 piece package | 534/S | 150X250mm conf. 1000 pz. 150X250mm 1000 pieces package |
| 408/S | rotoli 200 mm x 200 m conf. 2 pz. rolls 200 mm x 200 m 2 piece package | 536/S | 150X300mm conf. 1000 pz./ 150X300mm package 1000 pcs. |
| | | 538/S | 150X350mm conf. 1000 pz./ 150X350mm package 1000 pcs. |
| 410/S | rotoli 250 mm x 200 m conf. 2 pz. rolls 250 mm x 200 m 2 piece package | 540/S | 150X400mm conf. 1000 pz. 150X400mm 1000 pieces package |
| 412/S | rotoli 300 mm x 200 m conf. 2 pz. rolls 300 mm x 200 m 2 piece package | 550/S | 200X250mm conf. 1000 pz. 200X250mm 1000 pieces package |
| 414/S | rotoli 400 mm x 200 m conf. 1 pz. rolls 400 mm x 200 m 1 piece package | 540/S | 200X300mm conf. 1000 pz. 200X300mm 1000 pieces package |
| CODICE CODE | BUSTE AUTOSIGILLANTI AUTO-SEALING BAGS | | |
| | | 554/S | 200X350mm conf. 1000 pz./ 200X350mm package 1000 pcs. |
| | | 556/S | 200X400mm conf. 1000 pz./ 200X400mm package 1000 pcs. |
| 602/S | 90x230 mm (200 pz. per conf.) scatola da 12 conf. 90x230 mm (200 piece for package) 12 package box | 566/S | 250X350mm conf. 1000 pz./ 250X350mm package 1000 pcs. |
| | | 568/S | 250X400mm conf. 1000 pz./ 250X400mm package 1000 pcs. |
| 604/S | 140x250 mm (200 pz. per conf.) scatola da 8 conf. 140x250 mm (200 piece for package) 8 package box | 570/S | 250X450mm conf. 1000 pz./ 250X450mm package 1000 pcs. |
| | | 572/S | 300X500mm conf. 1000 pz./ 300X500mm package 1000 pcs. |
| 606/S | 190x330 mm (200 pz. per conf.) scatola da 6 conf. 190x330 mm (200 piece for package) 6 package box | 578/S | 300X400mm conf. 1000 pz./ 300X400mm package 1000 pcs. |
| | | 580/S | 300X500mm conf. 1000 pz./ 300X500mm package 1000 pcs. |
| 610/S | 300x450 mm (200 pz. per conf.) scatola da 2 conf. 300x450 mm (200 piece for package) 2 package box | 582/S | 300X550mm conf. 1000 pz./ 300X550mm package 1000 pcs. |
| | | 584/S | 400X500mm conf. 1000 pz./ 400X500mm package 1000 pcs. |
| | | 586/S | 400X600mm conf. 1000 pz. 400X600mm 1000 pieces package |
| CODICE CODE | BUSTE STERILIZZAZIONE PIATTE SALDATE SU 3 LATI FLAT STERILIZATION BAGS 3 WELDED SIDES | | |
| | | 588/S | 500X600mm conf. 1000 pz. 500X600mm 1000 pieces package |
| 514/S | 100X150mm conf. 1000 pz. 100X150mm 1000 pieces package | 590/S | 600X600mm conf. 1000 pz. 600X600mm 1000 pieces package |

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Testing the airtightness of the chamber (vacuum test)

Procedure: prepare the autoclave for work (preliminary warming up)-the chamber must be empty and closed start appropriate VT cycle (predisposed on the autoclave) wait for the end of the cycle and consult the printed report

Testing the steam-penetration through the indirect method (Bowie-Dick)

Procedure: prepare the Bowie Dick package by inserting the indicator place prepared package on the autoclave's central tray start standard sterilising cycle, by applying the following parameters: temperature: 134°C pressure: 2,1 bar extract and check indicator from the package at the end of the cycle if the indicator shows a homogenous change of colour from centre to border, the test can be considered to be correct. On the contrary, different and paler change of colour (or even if the colour doesn't change at all) indicate that the autoclave is not functioning correctly
Interventions: repeat the test if the result is negative if the second effort is positive, the autoclave is ready for use if the second test is still negative, the autoclave is not working correctly and thus it must be controlled and repaired by skilled technicians.

Testing the steam-penetration through the Helix Test method

Procedure: prepare the Helix test by inserting the indicator in the capsule at the edge of the tube put the test on the autoclave's central tray choose appropriate program and start-up sterilisation
Indications: check the capsule's indicator at the end of the process if the indicator shows a homogenous change of colour from centre to border, the test can be considered to be correct. On the contrary, different and paler change of colour, or even if the colour doesn't change at all (which can be usually noticed in the centre of the indicator), indicate that the autoclave is not functioning correctly because of the presence of air-bubbles

Chemical control

Procedure: During each process of sterilisation, it is necessary to check chemical indicators (whether they have changed colour or not) and eventual process-integrators that are used for all packages. On the contrary, it is necessary to establish the cause of the failed change of colour, prepare a new package and sterilise the instruments once again.
Integrators of process (applied to packages containing prosthetic materials) which changed colour must be saved and attached to

the section of chemical reports of biological control.

Biological control

It must be effected periodically.

Procedure: Insert the ampulla (which contains spores and agar agar) in the autoclave. Other materials can be inserted too. Choose appropriate cycle of sterilisation, then start up the process. Extract the ampulla from the chamber. Leave autoclave's doors open at least 15 minutes before taking it out. Handling the ampulla before it cools might break the glass and injure the operators. Thus, it is recommended the use of protective mask and gloves. It is necessary to keep an eye on the solution containing *Bacillus Sterothermophilus* 48 hours after delivery. This solution is used to control if the sterilising cycle has eliminated all spores within the indicator (negative result) or if they have survived (result of the test= *Bacillus Sterothermophilus*). It is necessary to repeat the test if, during the surveillance, one or more indicators detect the development of microorganism. In case of subsequent ascertainment of such development call skilled technician.

Interventions: If biological control is positive (i.e. presence of *Bacillus Sterothermophilus*), then:

- do not use your autoclave
- request intervention of service operators for special assistance and maintenance
- find (if possible) any other medical instrument which would give a negative result of the same process
- package once again the instruments mentioned above and repeat the cycle using another autoclave

After having repaired your autoclave, proceed with physical and microbiological verifications by means of vacuum test, Bowie Dick test, and biological checkups.

Autoclave can be put in operation only after having obtained negative results from all biological indicators that are used during the checkup-cycle and if physical parameters are all correct.

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
A correct sterilization, using the correct instruments, is useful in order to protect oneself from damages compensation claims.

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* The specifications given above are indicative. The operator must comply with regional or national laws/standards regulating sterilization.

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The tests certify, as a form of legal protection, that the equipment is valid and in working order and they must be filed and kept for 10 years, after which the professional's liability expires. Filing and keeping the tests carried out is a conscientious activity that may also be useful to the professional in the event of contentious procedure with patients.

As from today, filing is simple and inexpensive Tecno-Gaz S.p.A. has designed a new filing system called **STERILDOCS**. **STERILDOCS** is a system for filing all tests carried out on autoclaves:

- biological
- integrator
- bowie&dick
- helix test
- vacuum test

Sterildocs is a "nominal system" to which the serial number of the autoclave is attributed and is "fully equipped" for a year's use.

Filing is simple, rational and orderly and allows you to manage all the tests of each autoclave in an organized way.

Sterildocs is a product made by Tecno-Gaz, specialists in the sterilization field.

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
Operating protocols

"Guidelines", an extraordinary service that allows you to have written procedures.

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Dvd-video, a film that gives a practical demonstration of all the procedures and operating methodologies that operators must implement to obtain correct sterilization. Procedures for using the periodic tests (Bowie and Dick - Helix test - Vacuum test - biological tests - integrators).

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Aim of the process

Sterilization is the physical process that leads to the destruction of all forms of living microorganisms.

To ensure this result, specific physical conditions must be guaranteed that take into account the variability of the species of microorganisms potentially present on the device to be treated and above all their possible state:

- vegetative form
- sporogenous form

Spores are, in fact, by far the forms that are most resistant to sterilizing agents.

The best-known and most used sterilizing agent is heat, in particular wet heat in the form of steam.

If the steam is put under pressure, temperatures can be reached that are the sterilizing conditions of the penetrable materials and the surfaces exposed to the agent.

Steam is the safest, most economical, quickest and most harmless means of sterilization.

The devices that allow steam to be put under pressure are called autoclaves, or less properly sterilizers, equipped with a perfectly airtight chamber that can withstand high pressures.

Tecno-Gaz autoclaves

Tecno-Gaz has a wide range of autoclaves, able to satisfy all operating requirements and guarantee maximum safety, reliability and operating ergonomics.

For further information about these products, please contact Tecno-Gaz S.p.A. or an authorized dealer for the relevant commercial documents.

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All rights reserved. Variations can be done without notice. Tecno-Gaz S.p.A. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.

www.tecnogaz.com

 **TECNO-GAZ**
i n d u s t r i e s